

## The User's Manual—LMH SFM01

### Basic Information Introduction

#### Product Introduction

LMH SFM01 is serum-free media, without serum, protein, any animal-derived components and ect, formulated with high quality raw materials, high quality and stability, and no serum is required for the culture process, it suitable for high-density cell suspension expansion of LMH cell, supporting LMH rapid expansion, high-density culture and avian adenovirus vaccine production under a completely serum-free protein-free system.



#### Application Scope

LMH SFM01 medium is suitable for suspension domestication and large-scale serum-free suspension culture of LMH cells. It can be used in the development and production process of avian adenovirus vaccine and other biological products, and can support the development to large-scale production, but cannot be directly used in the human body or used as a medication.

#### Shipping, Storage and Validity Period

Product	Catalog No.	Storage	Shipping	Validity Period
LMH SFM01	LQ42, Liquid	2°C ~ 8°C, Protect from light	2°C ~ 8°C, Protect from light	6 months
LMH SFM01	DP42, Powder	2°C ~ 8°C, Protect from light	2°C ~ 8°C, Protect from light	24 months

#### Protocol for Hydration of Powder Medium

1. Fill the mixing container with purified water (20 ~ 30°C) at 100% of the final volume.
2. Slowly add 23.30 g/L of powder medium with gentle stirring. Mix for 20 ~ 30 minutes.
3. According to the dosage of  $0.24 \pm 0.01$  g/L, add 5 M NaOH solution to the above solution, and slowly adjust the pH to 6.0 ~ 6.5. Mix for 10 ~ 22 minutes, at this time the solution should be clear and transparent.
4. Slowly add  $2.00 \pm 0.05$  g/L  $\text{NaHCO}_3$  to the solution. Mix for 10 ~ 20 minutes.
5. If the pH value of the solution is not in the range of 7.0 ~ 7.4 at this time, adjust the pH to 7.0 ~ 7.4 with 5M NaOH

or 1M HCl solution.

6. Filter immediately the media with a 0.22 µm membrane filter.

### Quality Index of Powder and Liquid Media

Product Index	LMH SFM01 (LQ42), Liquid	LMH SFM01 (DP42), Powder
Appearance	Yellow, clear liquid	Light yellow or similar color powder
pH	7.0 ~ 7.6	7.0 ~ 7.4 (pre-filter)
Osmolality (mOsmol/kg)	280 ~ 340	280 ~ 340
Solubility	--	Dissolve well according to the protocol for hydration of powder medium
Endotoxin (EU/mL)	< 10	< 10
Sterility	Negative	--
Bioburden	--	Aerobic bacteria: < 200 CFU/g Molds and yeasts: < 200 CFU/g

## Reference Cell Culture Protocol

### Culture Conditions

	Parameter	Value
Culture volume	125 mL Shake flask	15 ~ 40 mL
	250 mL Shake flask	40 ~ 80 mL
	500 mL Shake flask	100 ~ 200 mL
	1000 mL Shake flask	200 ~ 300 mL
	Shake flask	25mm amplitude: 120 ~ 130 rpm
	Shake flask	50mm amplitude: 100 ~ 120 rpm
Culture environment	Seeding density	1.0×10 <sup>6</sup> cells/mL
	Incubation temperature	37°C
	Incubation CO <sub>2</sub> concentration	5%
	Incubation relative humidity	> 80% RH

## Cell Thawing

1. Pre-warm the medium (LMH SFM01) in 37°C water bath.
2. Spray the outside of the medium bottle with 75% alcohol and place the bottle into the bio-safety cabinet.
3. Thaw one vial at a time in 37°C water bath. Gently agitate the vial within 1 minutes until the ice in the vial melting.
4. Pipet the contents from the vial gently into a centrifuge tube containing 30 mL of pre-warmed medium (LMH SFM01).
5. Centrifuge 170 g ~ 190 g for 5 minutes. Discard the supernatant and re-suspend cells in 10 mL fresh pre-warmed medium (LMH SFM01), then adjust the cell density to  $1 \times 10^6$  cells/mL.
6. Sample 0.5 mL of cell suspension and analyze the viable cell density ( $\times 10^6$  cells/mL) and viability (%) of the sample using cell counter.
7. If the cell viability > 90%, incubate cells in the specified condition (refer to "culture conditions" table).

## Cell Passage

1. Pre-warm the medium (LMH SFM01) in 37°C water bath for 20 ~ 30 min.
2. Cells with viable cell density  $\geq 4 \times 10^6$  cells/mL, cell viability  $\geq 90\%$ , and in the middle of logarithmic growth phase were selected for passage.
3. According to the seed cell density of  $1.0 \times 10^6$  cells/mL, calculate the amount of total number of seed cells.
4. Seed cells at  $1.0 \times 10^6$  cells/mL in the flask and add a certain volume of pre-warmed fresh medium.
5. Incubate cells in the specified environment condition (refer to "culture conditions" table).
6. Passage cells with fresh LMH SFM01 medium according to the above steps once every  $48 \pm 3$  hours.
7. If the viable cell density is less than  $4.0 \times 10^6$  cells/mL or the viability is lower than 90% before passaging, the cells need to be centrifuged at 170 g ~ 190 g (about 800 rpm ~ 1000rpm) for 5 minutes. Carefully remove the spent media, then resuspend cells with preheated LMH SFM01 medium, passage cells after sampling and counting.

## Adaptation

### Direct Adaptation

1. For cells can direct adapt, transfer cells suspension cultures into LMH SFM01 directly, and the seed cell density refer to the cell passage procedure.
2. Cell passage until cell expression steadily.

3. When VCD reaches  $4 \times 10^6$  cells/mL and  $> 90\%$  viability ( $48 \pm 3$  hours) . At this point, the cells had been successful adapted.

### **Sequential Adaptation**

For cells growing in 5 ~ 10% serum or SFM media. Sequential adaptation should be performed.

1. Seed cells at  $1.0 \times 10^6$  cells/mL in original cell culture media.
2. Sample and cell count every day until the VCD reaches  $4 \times 10^6$  cells/mL.
3. Seed cell density at  $1.0 \times 10^6$  cells/mL, subculture cells into stepwise increasing ratios of complete LMH SFM01 medium to original medium with each subsequent passage (25:75, 50:50, 75:25, 90:10, 100:0).
4. When VCD reaches  $4 \times 10^6$  cells/mL and  $\geq 90\%$  viability ( $48 \pm 3$  hours) . At this point, the cells had fully adapted to LMH SFM01 media.

### **Cell Cryopreservation**

1. Prepare cells, harvesting in mid-log phase of growth with viability  $> 90\%$ .
2. Sample and cell counting, calculate the required volume of cell freeze solution to give a final density of  $2.5 \sim 3.5 \times 10^7$  cells/mL.
3. Prepare the cell freeze solution: 93% LMH SFM01 + 7% DMSO, store at  $4^\circ\text{C}$ .
4. Centrifuge 170 g ~ 190 g for 5 minutes, discard the supernatant and re-suspend cells with the cell freeze solution.
5. Immediately dispense aliquots of cells suspension into cryovials according to the specific needs of the project.
6. Achieve cryopreservation in an automated or manual controlled rate freezing apparatus following standard procedures ( $1^\circ\text{C}$  decrease per minute).
7. Transfer to liquid nitrogen tank for storage.