

BHK Serum-free Medium

Product Introduction

BHK SFM01 is serum-free, protein free, animal origin-free medium. It is suitable for BHK cell suspension domestication and scale up serum-free suspension culture without the addition of serum during culture. It supports rapid BHK expansion, high density culture and virus vaccine production in a completely serum-free and protein-free system. It can be used in the R&D and production process of pseudo rabies virus, etc., and can support from R&D to scale up production.



Product Information

Product	Catalogue	Characteristic	Suitable Cell Line & Related Derivatives
BHK SFM01 Serum-free Medium	LQ16-1000 (Liquid)	serum-free, protein-free and animal origin free	BHK, BHK21
	DP16,50L/100L (Dry Powder)		

Product Advantage

- Serum-free, protein-free and animal origin free;
- (BHK Cells)Rapid adaptation to suspension;
- Support BHK cells rapid proliferation and high density culture;
- Support efficient expansion of rabies virus;
- Suitable for large scale bioreactor culture processes, supporting batch, fed-batch and perfusion modes;
- Highly quality raw materials to formulate the product with high quality and stability;
- cGMP standard production conditions, using continuous needle milling equipment for power medium production to ensure good product consistency of different batches.

