

产品规格 / PRODUCT SPECIFICATIONS
---Vero SFM01 Pro

产品名称 Product Name	Vero SFM01 Pro With Glucose (3.90g/L), L-Glutamine, Phenol red
产品编号 Product Code	DP22
产品特点 Characteristics	<ul style="list-style-type: none"> ● 化学组分限定、无动物源培养基，适合VERO贴壁静态培养及微载体悬浮规模化培养，适用于疫苗大规模生产表达。 ● Chemical-defined, animal origin-free medium optimized for Vero adherent culture and large-scale suspension cultures in microcarrier, suitable for vaccine large-scale manufacturing applications.
形式 Format	粉末，非无菌 / Powder, Non-sterile
包装规格 Pack Size	50L/100L (16.85 g/L)
存放条件 Storage	2°C ~ 8°C，避光 / 2°C ~ 8°C, Protected from light
有效期 Shelf life	24 个月 / 24 Months