

产品规格 / PRODUCT SPECIFICATIONS
---Tobitec FA Pro

产品名称 Product Name	Tobitec FA Pro With Glucose (75.00g/L) Without L-Glutamine
产品编号 Product Code	DPF02-A
产品特点 Characteristics	<ul style="list-style-type: none"> ● 化学成分确定的补料培养基，无动物源性成分，无蛋白。提高CHO细胞生长和生产性能专门设计的一种高效补料配方，可以满足细胞高密度培养和产物高表达时的营养需求。结合细胞株特性，采用合理的补料策略可以明显改善细胞的生长和抗体的表达量以及质量。 ● Chemical defined, animal origin-free, protein-free medium. A highly effective feeding formulation specially designed to improve the growth and production performance of CHO cells can meet the nutritional needs of cell high-density culture and high expression of products. Combined with the characteristics of cell lines, a reasonable feeding strategy can significantly improve cell growth and antibody expression and quality.
形式 Format	粉末，非无菌 / Powder, Non-sterile
包装规格 Pack Size	5L/10L (181.04 g/L)
存放条件 Storage	2°C ~ 8°C，避光 / 2°C ~ 8°C, Protected from light
有效期 Shelf life	24 个月 / 24 Months