

产品规格 / **PRODUCT SPECIFICATIONS**

---**Tobitec FB**

<p>产品名称 <b>Product Name</b></p>	<p>Tobitec FB Without Glucose, L-Glutamine</p>
<p>产品编号 <b>Product Code</b></p>	<p>DPF01-B</p>
<p>产品特点 <b>Characteristics</b></p>	<ul style="list-style-type: none"> <li>● 化学成分确定的补料培养基，无动物源性成分，无蛋白。提高CHO细胞生长和生产性能专门设计的一种高效补料配方，可以满足细胞高密度培养和产物高表达时的营养需求。结合细胞株特性，采用合理的补料策略可以明显改善细胞的生长和抗体的表达量以及质量。</li> <li>● Chemical defined, animal origin-free, protein-free medium. A highly effective feeding formulation specially designed to improve the growth and production performance of CHO cells can meet the nutritional needs of cell high-density culture and high expression of products. Combined with the characteristics of cell lines, a reasonable feeding strategy can significantly improve cell growth and antibody expression and quality.</li> </ul>
<p>形式 <b>Format</b></p>	<p>粉末，非无菌 / Powder, Non-sterile</p>
<p>包装规格 <b>Pack Size</b></p>	<p>10L/20L (94.50 g/L)</p>
<p>存放条件 <b>Storage</b></p>	<p>2°C ~ 8°C，避光 / 2°C ~ 8°C, Protected from light</p>
<p>有效期 <b>Shelf life</b></p>	<p>24 个月 / 24 Months</p>