

产品规格 / PRODUCT SPECIFICATIONS**--- BHK SFM01**

产品名称 Product Name	BHK SFM01 With Glucose (8.70g/L), L-Glutamine
产品编号 Product Code	DP16
产品特点 Characteristics	<ul style="list-style-type: none">● 无动物源、无血清培养基，适合BHK21细胞无血清生长培养。● Animal origin-free, serum-free medium suitable for the culture of BHK21 cell.
形式 Format	粉末，非无菌 / Powder, Non-sterile
包装规格 Pack Size	50L/100L (24.724 g/L)
存放条件 Storage	2°C ~ 8°C，避光 / 2°C ~ 8°C, Protected from light
有效期 Shelf life	24 个月 / 24 Months