

产品规格 / PRODUCT SPECIFICATIONS
---Insect SFM01 Pro

产品名称 Product Name	Insect SFM01 Pro With Glucose (11.00g/L), L-Glutamine
产品编号 Product Code	DP20
产品特点 Characteristics	<ul style="list-style-type: none"> ● 无血清,无动物源组分, 专门为Insect细胞系开发, 可用于Insect悬浮高密度培养即杆状病毒转染表达, 适用于大规模生产应用。 ● 联合 Insect TE02使用。 ● No serum, no animal origin components, specially developed for Insect cell line, can be used for suspension high density culture of Insect namely baculovirus transfection expression, suitable for large-scale production application. ● Using combined with Insect TE02.
形式 Format	粉末, 非无菌 / Powder, Non-sterile
包装规格 Pack Size	50L (43.40 g/L)
存放条件 Storage	2°C ~ 8°C, 避光 / 2°C ~ 8°C, Protected from light
有效期 Shelf life	24 个月 / 24 Months