

## 产品规格 / PRODUCT SPECIFICATIONS

### ---T淋巴细胞无血清培养基 T Lymphocyte Serum-Free Culture Medium

<b>产品名称</b> <b>Product Name</b>	T淋巴细胞无血清培养基 T Lymphocyte Serum-Free Culture Medium
<b>产品编号</b> <b>Product Code</b>	CM012-500/1000
<b>产品特点</b> <b>Characteristics</b>	<ul style="list-style-type: none"><li>● 无血清，无动物源组分，适用于T细胞（T、CAR-T及TCR-T等）体外无血清扩增培养。</li><li>● Serum-Free, animal origin-free, suitable for serum-free expansion culture of T cells (T, CAR-T, TCR-T, etc.) in vitro.</li></ul>
<b>形式</b> <b>Format</b>	液体，无菌 / Liquid, Asepsis
<b>包装规格</b> <b>Pack Size</b>	500, 1000mL/瓶
<b>存放条件</b> <b>Storage</b>	2°C ~ 8°C，避光 / 2°C ~ 8°C, Protected from light
<b>有效期</b> <b>Shelf Life</b>	12 个月/12 Months