

产品规格 / PRODUCT SPECIFICATIONS
---HEK293 CDM26 Pro

产品名称 Product Name	HEK293 CDM26 Pro With Glucose (7.00g/L), L-Glutamine
产品编号 Product Code	DP26-2
产品特点 Characteristics	<ul style="list-style-type: none"> ● 化学成分限定,适应于HEK293 (293F、293T、293H、Expi293等) 悬浮高密度培养, 蛋白及慢病毒等转染表达, 适合大规模生产应用。 ● Chemical defined media, animal origin-free medium specially designed for HEK293.It can be used for the growth of HEK293 (293F,293T,293H,Expi293) suspension cells high-density culture, adenovirus and proteins expression, suitable for large-scale manufacturing applications
形式 Format	粉末, 非无菌 / Powder, Non-sterile
包装规格 Pack Size	50L/100L (23.00 g/L)
存放条件 Storage	2°C ~ 8°C, 避光 / 2°C ~ 8°C, Protected from light
有效期 Shelf life	24 个月 / 24 Months