

产品规格 / PRODUCT SPECIFICATIONS
---CHO CDM34

产品名称 Product Name	CHO CDM34 With Glucose (6.00g/L), Sodium pyruvate Without L-Glutamine, HEPES
产品编号 Product Code	DP34
产品特点 Characteristics	<ul style="list-style-type: none"> ● 化学成分限定，无动物源组分，无蛋白，适合CHO系列细胞悬浮高密度大规模培养，可提高蛋白及抗体的表达与质量。 ● Chemical defined, animal origin-free, protein-free medium optimized for the growth of high-density and large-scale suspension cultures of CHO series cells, and can improve the expression and quality of antibodies and proteins.
形式 Format	粉末，非无菌 / Powder, Non-sterile
包装规格 Pack Size	50L/100L (21.90 g/L)
存放条件 Storage	2°C ~ 8°C，避光 / 2°C ~ 8°C, Protected from light
有效期 Shelf life	24 个月 / 24 Months