

产品规格 / **PRODUCT SPECIFICATIONS**

---**Insect SFM01 Pro**

产品名称 Product Name	Insect SFM01 Pro With Glucose (11.00g/L), L-Glutamine
产品编号 Product Code	DP20
产品特点 Characteristics	● 无血清,无动物源组分，专门为Insect细胞系开发，可用于Insect悬浮高密度培养即杆状病毒转染表达，适用于大规模生产应用。 ● 联合 Insect TE02使用。 ● No serum, no animal origin components, specially developed for Insect cell line, can be used for suspension high density culture of Insect namely baculovirus transfection expression, suitable for large-scale production application. ● Using combined with Insect TE02.
形式 Format	粉末，非无菌 / Powder, Non-sterile
包装规格 Pack Size	50L (43.40 g/L)
存放条件 Storage	2℃ ~ 8℃，避光 / 2℃ ~ 8℃, Protected from light
有效期 Shelf life	36 个月 / 36 Months