

**产品规格 / PRODUCT SPECIFICATIONS**
**--- MEM $\alpha$ , nucleosides**

<b>产品名称</b> <b>Product Name</b>	MEM $\alpha$ , nucleosides With Glucose (3.15g/L), L-Glutamine, Sodium pyruvate Without NaHCO <sub>3</sub> , HEPES, Phenol red
<b>产品编号</b> <b>Product Code</b>	DP36A
<b>产品特点</b> <b>Characteristics</b>	<ul style="list-style-type: none"> <li>● 支持多种哺乳动物细胞系的生长。</li> <li>● Basal medium for supporting the growth of many different mammalian cell lines.</li> </ul>
<b>形式</b> <b>Format</b>	粉末，非无菌 / Powder, Non-sterile
<b>包装规格</b> <b>Pack Size</b>	50L/100L (10.12 g/L)
<b>存放条件</b> <b>Storage</b>	2°C ~ 8°C, 避光 / 2°C ~ 8°C, Protected from light
<b>有效期</b> <b>Shelf life</b>	36 个月 / 36 Months